## Antibiotic Antimycotic Solution (100X)

Contains 10,000 units/mL of penicillin  
10 mg/mL streptomycin sulfate  
25 µg/mL amphotericin B  
in 0.85% sodium chloride

For gram-positive, gram-negative bacteria, yeasts, and molds

Catalog Number **LS 203-01**

Storage Temperature -5~-20°C

### Product Description

Microbial contamination results in the lethal errors in cell and tissue culture. Animal cell divides into two daughter cells in several hours to days, however microbes such as bacteria, fungi, or yeasts are doubled in decades of minute to hours. Hence, overgrowing of these microbes inhibits the growth, proliferation, and function of animal cells. Microbial components or products can result the toxic effect to animal cells. Generally, in culture medium, antibiotics are added for inhibition of microbial growth. Though there are many types of antibiotics used in cell culture media, penicillin G, streptomycin, and amphotericin B are widely used.

**LS 203-01** may be used as the antibiotic and antimycotic solution for preventing the contamination of bacteria, yeasts, and molds. The generally used concentrations of these antibiotic and antimycotic materials are 100,000 units/L for penicillin, 100 mg/L for streptomycin, and 0.25 mg/L for amphotericin B. Each concentration of antibiotic and antimycotic materials can be obtained by 1/100 dilution of this solution in medium.

### Biological Performance Characteristics

The growth-promoting capacity of the antibiotics is tested in a medium containing 10% FBS using an appropriate cell line(s). Growth rates are examined through three subculture generations and compared with parallel cultures grown in standardized control medium. Cells are counted and growth is plotted as a logarithmic function of time in culture, and seeding efficiencies, doubling time, and final cell densities are determined. During the testing period, cultures are examined microscopically for a typical morphology and evidence of cytotoxicity.

The concentration of the antibiotics should be reduced to 1/2~1/10 for cell culture using serum free media.

### Precautions

For *In Vitro* Use Only

Product Profile	
Appearance	Clear colorless solution
pH	6.1 ~ 6.7
Osmolality	299 ~ 347 mOsm/kg H <sub>2</sub> O
Endotoxin	≤ 10 EU/ml
Sterility	Sterilized by 0.2 µm filtration system. Sterility tests are performed in accordance with protocols described in USP.

### References

Perlman, D., Use of antibiotics in cell culture media, in *Methods in Enzymology*, Jakoby, W. and Pastan, I. H., eds, Academic Press, New York, NY, 1979, Vol. LVIII, 112.  
Data for *Biochemical Research*, 3<sup>rd</sup> Edition, Dawson, R. M. C. et al., eds., Oxford University Press, Inc. New York, 1986, p 297.  
Alberts, B. et al., Basic genetic mechanisms, in *Molecular Biology of the Cell*, 3<sup>rd</sup> Edition, Garland Publishing, Inc., New York, NY, 1994, p 240.  
Prophylactic use of antibiotics in cells and tissues with a high risk of microbial contamination, in *Cell and Tissue Culture: Laboratory Procedures*, Doyle, A. et al., eds., John Wiley and Sons, England, 1996, Vol. 1, p 2A:4.2.

Antibiotics	Target Microbes*	Actions
Penicillin	G (+)	Interferes with the final stage of synthesis of the bacterial cell wall
Streptomycin	G (+, -)	Binds to 30S subunit to cause misreading
Amphotericin B	Y, M	Interferes with the permeability of cell membrane of sensitive fungi by binding sterols

\* G (+), Gram-positive bacteria  
G (+, -), Gram-positive and Gram-negative bacteria  
Y, Yeasts  
M, Molds

### Storage/Stability

Antibiotic antimycotic solution should be stored at -5~-20°C in the dark and can be stable for 3 days during cultivation of cells at 37°C. After 3 days, medium should be refreshed to a new medium. Deterioration of the liquid may be recognized by (1) precipitate or particulate matter throughout the solution, (2) cloudy appearance, (3) color change, and/or (4) pH change. The nature of supplements added may affect storage conditions and shelf life of the medium. Product label bears expiration date.

## Antibiotic Antimycotic Solution (100X)

Contains 10,000 units/mL of penicillin  
10 mg/mL streptomycin sulfate  
25 µg/mL amphotericin B  
in 0.85% sodium chloride

For gram-positive, gram-negative bacteria, yeasts, and molds

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For *In Vitro* Use Only

Product Profile	
Appearance	Clear colorless solution
pH	6.1 ~ 6.7
Osmolality	299 ~ 347 mOsm/kg H <sub>2</sub> O
Endotoxin	≤ 10 EU/ml
Sterility	Sterilized by 0.2 µm filtration system. Sterility tests are performed in accordance with protocols described in USP.

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Perlman, D., Use of antibiotics in cell culture media, in *Methods in Enzymology*, Jakoby, W. and Pastan, I. H., eds, Academic Press, New York, NY, 1979, Vol. LVIII, 112.  
Data for Biochemical Research, 3<sup>rd</sup> Edition, Dawson, R. M. C. et al., eds., Oxford University Press, Inc. New York, 1986, p 297.  
Alberts, B. et al., Basic genetic mechanisms, in *Molecular Biology of the Cell*, 3<sup>rd</sup> Edition, Garland Publishing, Inc., New York, NY, 1994, p 240.  
Prophylactic use of antibiotics in cells and tissues with a high risk of microbial contamination, in *Cell and Tissue Culture: Laboratory Procedures*, Doyle, A. et al., eds., John Wiley and Sons, England, 1996, Vol. 1, p 2A:4.2.

Antibiotics	Target Microbes*	Actions
Penicillin	G (+)	Interferes with the final stage of synthesis of the bacterial cell wall
Streptomycin	G (+, -)	Binds to 30S subunit to cause misreading
Amphotericin B	Y, M	Interferes with the permeability of cell membrane of sensitive fungi by binding sterols

\* G (+), Gram-positive bacteria  
G (+, -), Gram-positive and Gram-negative bacteria  
Y, Yeasts  
M, Molds

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